


## AvaLight-DHc Compact Deuterium Halogen Light Source

The AvaLight-DHc is a combined deuterium-halogen light source, to be used for UV/VIS/NIR applications. The output energy of the AvaLight-DHc is relatively low. The use is therefore recommended in transmission setup with large diameter fibers. The light source emits light from 200 to 2500 nm and has an SMA connector to easily couple into our range of fiber optics. The AvaLight-DHc has an integrated TTL shutter, that can be used for auto-save dark/ lamp off in combination with AvaSoft (extra IC-DB15-2 or IC-DB26-2 needed).

### AvaLight-DHc



Figure 8 Spectral output of AvaLight-DHc



Optionally, the AvaLight-DHc can be delivered in Rackmounted version, to be fully integrated in the 19" rackmount or 9.5" desktop housing.

A direct attach cuvette holder CUV-DHc/XE/LED (see section accessories) is available for fluorescence or absorbance measurements.

### Technical Data


	Deuterium Light Source	Halogen Light Source
<b>Wavelength Range</b>	200 - 400 nm	400 - 2500 nm
<b>Stability</b>	< 1 mAU	< 1 mAU
<b>Warm-up time</b>	8 min	1 min
<b>Drift</b>	< 0.25% / h	< 0.25% / h
<b>Optical Power in 600 <math>\mu</math>m fiber</b>	0.2 $\mu$ Watt	7 $\mu$ Watt
<b>Lamp Lifetime</b>	1000 hours	2000 hours
<b>Temp. Range</b>	5°C - 35°C	
<b>Power Supply</b>	12VDC / 450mA	
<b>Dimensions</b>	175 x 110 x 44 mm	

### ORDERING INFORMATION

<b>AvaLight-DHc</b>	Compact Deuterium Halogen Light Source with TTL Shutter
<b>IC-DB15-2</b>	Interface cable AvaSpec-USB1 platform to AvaLight-DHc-TTL-shutter
<b>IC-DB26-2</b>	Interface cable AvaSpec-USB2 platform to AvaLight-DHc-TTL-shutter
<b>AvaLight-DHc-RM</b>	Rackmount Compact Deuterium Halogen Light Source with TTL Shutter
<b>AvaLight-DHc-B</b>	Compact Deuterium Halogen Replacement Bulb
<b>CUV-DHc/XE/LED</b>	Direct attach cuvette holder for AvaLight-DHc/XE/LED
<b>PS-12V/1.0A</b>	Power supply 100-240VAC/12VDC, 1.0A for AvaLight-DHc

## AvaLight-HAL Tungsten Halogen Light Source

### AvaLight-HAL



The AvaLight-HAL is a compact stabilized halogen fan-cooled light source that can be used for the visible and the near infrared range. The AvaLight-HAL features adjustable focusing of the SMA connector to maximize light coupling into a fiber or fiber bundle with a diameter of up to 600  $\mu\text{m}$ . A filter slot accepts 1" round or 2" x 2" square filters up to 3 mm thick. The lamp stability is achieved by a current stabilization. A fan regulates the airflow around the heatsink to optimize the operation temperature. The AvaLight-HAL needs an extra PS-24V/1.25A 24 VDC power adapter.

The SMA-connector input into any fiber can be optimized by changing the focus. Bulb replacement is easy. With an internal jumper the optical output energy can be controlled. At "low" setting the source acts as a long life source with over 4000 hrs life time. At "medium" setting the color temperature goes up and the expected life time is about 2000 hrs. The "high" setting gives max output in the blue range, but reduces bulb life time to ca. 1000 hrs (see fig 9).


The AvaLight-HAL-S has an internal TTL shutter, that can be controlled from the AvaSpec, so the auto-save dark option in the AvaSoft software can be used (extra IC-DB15-2 or IC-DB26-2 needed).

The filter holder can be easily replaced by a direct attach cuvette holder CUV-HAL (see section accessories) useful for fluorescence or absorbance measurements.

Optionally, the AvaLight-HAL(-S) can be delivered in Rackmounted version, to be fully integrated in the 19" rack-mount or 9.5" desktop housing.

New is the AvaLight-HAL-S-IND halogen light source with a heavy duty industrial power connector and 24V power adapter included. We also offer to implement the new heavy duty industrial connector in your existing AvaLight-HAL-S and power supply as an upgrade.

Figure 9 Spectral output of AvaLight-HAL



## Technical Data

	AvaLight-HAL (standard)	AvaLight-HAL (long life)	AvaLight-HAL (high power)
<b>Wavelength Range</b>	360-2500 nm		
<b>Stability</b>	± 0.1%		
<b>Time to stabilize</b>	Ca. 15 min.		
<b>Output to bulb</b>	12.7 VDC/ 0.9A	11.3 VDC/ 0.8A	14.1 VDC/ 1.0A
<b>Bulb Life</b>	2000 hrs	> 4000 hrs	< 1000 hrs
<b>Optical power 200µm fiber</b>	0.5 mWatt	0.35 mWatt	0.7 mWatt
<b>Optical power 600µm fiber</b>	4.5 mWatt	3.2 mWatt	6 mWatt
<b>Optical power 1000µm fiber</b>	10 mWatt	7 mWatt	14 mWatt
<b>Bulb Color Temperature</b>	2.850 K	2.700 K	3.000 K
<b>Power requirement</b>	24 VDC / 1.25A		
<b>Temperature range</b>	0-70 °C		
<b>Dimensions, weight</b>	132 x 110 x 44 mm, 490 grams		

### Separate 50x50mm filters to install in AvaLight-HAL (-S)

<b>GL-WG305-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 305 nm
<b>GL-KG3-3</b>	Separate 50 x 50 x 3 mm bandpass filter, transparent > 325 nm and < 700 nm
<b>GL-BG28-3</b>	Separate 50 x 50 x 3 mm bandpass filter, transparent > 360 nm and < 500 nm
<b>GL-GG385-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 385 nm
<b>GL-GG475-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 475 nm
<b>GL-OG515-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 515 nm
<b>GL-OG550-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 550 nm
<b>GL-OG590-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 590 nm
<b>GL-NG9-1</b>	Separate 50 x 50 x 1 mm Neutral Density filter (transmission ca. 10% 400-1100nm)
<b>GL-NG9-2</b>	Separate 50 x 50 x 2 mm Neutral Density filter (transmission ca. 1% 400-1100nm)
<b>GL-NG9-3</b>	Separate 50 x 50 x 3 mm Neutral Density filter (transmission ca. 0.1% 400-1100nm)

**More filter types available, please contact us for ordering information**

### ORDERING INFORMATION

<b>AvaLight-HAL</b>	10W Tungsten Halogen Lamp, fan-cooled, needs extra PS-24V/1.25A power supply
<b>AvaLight-HAL-S</b>	10W Tungsten Halogen Lamp, fan-cooled, incl. TTL shutter, needs extra PS-24V/1.25A power supply
<b>AvaLight-HAL-S-IND</b>	10W Tungsten Halogen Lamp, fan-cooled, incl. TTL shutter, including industrial connector and PS-24V/1.25A power supply
<b>AvaLight-HAL-S-RM</b>	Rack mounted version of the 10W Tungsten Halogen Lamp, fan-cooled, incl. TTL shutter
<b>IC-DB15-2</b>	Interface cable AvaSpec-USB1 platform to AvaLight-HAL-S
<b>IC-DB26-2</b>	Interface cable AvaSpec-USB2 platform to AvaLight-HAL-S
<b>AvaLight-HAL-B</b>	10W Tungsten Halogen Replacement bulb for AvaLight-HAL, AvaLight-HAL-S
<b>HAL-S-IND-UPGRADE</b>	Upgrade from HAL-S to HAL-S-IND, incl. industrial connector please return with PS-24V/1.25A
<b>CUV-HAL</b>	Direct attach cuvette holder for AvaLight-HAL(-S)
<b>PS-24V/1.25A</b>	Power supply 100-240VAC/24VDC, 1.25A, necessary for AvaLight-HAL

## AvaLight-D(H)-S Deuterium-Halogen Light Sources

### AvaLight-DH-S



The AvaLight-DH-S is a combined Deuterium and Halogen light source, which can be used for UV/VIS/NIR applications. The AvaLight-DH-S has an SMA905 connector to be used with fiber optic cables and bundles. For optimal coupling an adjustable focus lens assembly is included. The light source supplies a continuous spectrum with high efficiency and highest stability in the UV, visible and near-infrared range from 215-2500 nm.

The AvaLight-DH-S was developed based on the shine-through principle, in which the halogen light is focused through a small diameter aperture in the deuterium bulb. The AvaLight-DH-S-BAL has a dichroic beam splitter to give a much more balanced spectrum from 200-2500nm.

For UV only applications the AvaLight-D-S (only Deuterium) source is available.

The standard AvaLight-D(H)-S light sources have a special UV longpass filter (>220nm) implemented to protect the fibers from solarizing. For applications that require a spectral range below 220 nm a Deep UV bulb is available, in which no UV longpass filter is included. For these Deep-UV measurements (from 190nm) the AvaLight-D(H)-S-DUV may be ordered. New is a long-life and more stable version of the Deuterium lamp (AvaLight-D-B-DUV-LL).

For all Deuterium Light Sources solarization resistant fibers (-SR) are recommended, for the DUV lamp special solarization resistant fibers (-SR) are required (see the fiber optic section).

The output of the AvaLight-DH-S is optimized for fibers (bundles) with a diameter up to 600  $\mu\text{m}$ , for larger diameter fiber bundles the AvaLight-DH-S SMA connector focal point can be manually adapted to uniformly fill a larger core diameter fiber cable or bundles.

In all AvaLight-Deuterium sources an integrated TTL-shutter and filter holder for filters of up to 50x50x6.5 mm are implemented.

The AvaLight-DH-S internal TTL shutter can be controlled from the AvaSpec, so the auto-save dark option in the AvaSoft-FULL software can be used (extra IC-DB15-2 or IC-DB26-2 needed).

The filter holder can be easily replaced by a direct attach cuvette holder CUV-DHS (see section accessories) useful for fluorescence or absorbance measurements.

Figure 10 Spectral output AvaLight-D-S





Figure 11 Spectral output AvaLight-DH-S



## Technical Data

	Deep UV Deuterium	Deuterium (Standard)	Deuterium (long-life DUV)	Halogen Lamp
<b>Wavelength Range</b>	190-400nm	215-400nm	185-400nm	360-2500nm
<b>Warm-up Time</b>	30 min.	30 min.	30 min.	20 min.
<b>Lamp Power</b>	78W / 0.75A	78W / 0.75A	78W / 0.75A	5W / 0.5A
<b>Lamp Lifetime</b>	1000 hrs	1000 hrs	2000 hrs	1000 hrs
<b>Noise (AU)</b>	$2 \times 10^{-4}$	$2 \times 10^{-4}$	$2 \times 10^{-5}$	$10^{-4}$
<b>Max. drift</b>	$\pm 0.5\%/hrs$	$\pm 0.5\%/hrs$	$\pm 0.5\%/hrs$	$\pm 0.1\%/hrs$
<b>Color Temperature</b>	-	-	-	3000 K
<b>Optical Power in 200<math>\mu</math>m fiber</b>	11 $\mu$ W	7 $\mu$ W	11 $\mu$ W	43 $\mu$ W
<b>Optical Power in 600<math>\mu</math>m fiber</b>	72 $\mu$ W	61 $\mu$ W	72 $\mu$ W	239 $\mu$ W
<b>Optical Power in 1000<math>\mu</math>m fiber</b>	206 $\mu$ W	166 $\mu$ W	206 $\mu$ W	354 $\mu$ W
<b>Power consumption</b>	90 Watt (190Watt for heating D-Lamp 4-5 sec.)			
<b>Power Requirements</b>	100-240VAC 50/60 Hz			
<b>Dimensions / Weight</b>	315 x 165 x 140 mm / ca 5 kg.			

For a table of separate 50x50mm filters to install in AvaLight-D(H)-S see AvaLight-HAL.

### ORDERING INFORMATION

<b>AvaLight-D-S</b>	Deuterium light source, 215-400 nm, incl. TTL shutter, -SR fibers recommended
<b>AvaLight-DH-S</b>	Deuterium-Halogen light source, 215-2500 nm, incl. TTL shutter, -SR fibers recommended
<b>AvaLight-D-S-DUV</b>	Deep UV deuterium light source, 190-400 nm, incl. TTL shutter, -SR fibers needed
<b>AvaLight-DH-S-DUV</b>	Deep UV deuterium-halogen light source, 190-2500 nm, incl. TTL shutter, -SR fibers needed
<b>AvaLight-DH-S-DUV-LL</b>	Deep UV deuterium-halogen light source, long-life 2000 hrs, 190-2500 nm, incl. TTL shutter, -SR fibers needed!
<b>IC-DB15-2</b>	Interface cable AvaSpec-USB1 platform to AvaLight-D(H)S
<b>IC-DB26-2</b>	Interface cable AvaSpec-USB2 platform to AvaLight-D(H)S
<b>AvaLight-D-B</b>	Replacement deuterium bulb for AvaLight-D/AvaLight DH light source
<b>AvaLight-D-B-DUV</b>	Replacement deep UV deuterium bulb for AvaLight-D(H)-S-DUV light source
<b>AvaLight-D-B-DUV-LL</b>	Long-life 2000 hrs. Replacement deep UV deuterium bulb for AvaLight-D/AvaLight DH light source
<b>AvaLight-DH-B</b>	Replacement halogen bulb for AvaLight-DH light source
<b>CUV-DHS</b>	Direct attach cuvette holder for AvaLight-D(H)S