

Avantes Deuterium light sources are known for their stable output and are used for UV absorption or reflection measurements. These can also be used as irradiance calibration sources due to their high-stability. The standard AvaLight-DH-S mixes the Halogen light with the Deuterium light, thus producing a wide spectral range light source. The output spectrum of Deuterium light sources exhibits several peaks, with a prominent peak at 656 nm. The AvaLight-DH-S-BAL incorporates a dichroic beam splitter installed to minimize these peaks, providing a smooth spectrum from 200-2500 nm.

Our pulsed Xenon light source is used in applications where a long lifetime and high output power is needed, such as in fluorescence measurements. This is an affordable UV source, but the spectral output is not as smooth and continuous as the AvaLight Halogen and Deuterium light sources. LED light sources provide high power at a precise wavelength. A typical application for AvaLight-LED sources is fluorescence. They provide long lifetime, short warm-up time and high-stability.

For wavelength calibration Avantes offers a variety of sources including Argon, Mercury-Argon, Neon, Zinc and Cadmium. All Avantes spectrometers are factory wavelength calibrated and do not require recalibration as they have fixed slits and optics. For those customers who wish to do their own calibrations, the AvaLight-CAL light sources can be used for recalibration purposes. For auto-calibration AvaSoft-Full provides a calibration procedure to make this easy.

New to the AvaLight illumination sources family is the AvaLight-LDXE, a very high brightness light source. It has an extremely long lamp life, a spectral range from 200 nm to 1100 nm and a high-stability.

**Table 8 Overview light sources**

	Wavelength Range	Type	Principle	Product
<b>Color / VIS / NIR</b>	360-2500 nm	Tungsten Halogen	Continuous	AvaLight-HAL(-S)-MINI
<b>DUV</b>	190-400 nm	Deuterium	Continuous	AvaLight-D-S-DUV
<b>UV</b>	215-400 nm	Deuterium	Continuous	AvaLight-D-S
<b>UV/VIS/NIR refl./abs</b>	215-2500 nm	Deuterium/Halogen	Continuous	AvaLight-DH-S-(BAL)
<b>UV/VIS/NIR absorption</b>	200-2500 nm	Deuterium/Halogen	Continuous	AvaLight-DHc
<b>UV/VIS</b>	200-1000 nm	Xenon	Pulsed	AvaLight-XE
<b>Fluorescence</b>	Multiple possible	LED	Continuous	AvaLight-LED
<b>Wavelength Calibration</b>	253-1704 nm	Mercury-Argon Neon / Argon	Continuous	AvaLight-CAL
	200-700 nm	Zinc/Cadmium	Continuous	AvaLight-CAL-CAD/Zinc
<b>Irradiance Calibration</b>	360-2500 nm	Tungsten Halogen	Continuous	AvaLight-HAL-CAL-MINI
	200-1100 nm	Deuterium/Halogen	Continuous	AvaLight-DH-(BAL)-CAL
<b>Radiance Calibration</b>	360-2500 nm	Tungsten Halogen	Continuous	AvaSphere-50-LS-HAL-CAL




Figure 7 spectral distribution of different Light Sources

# AvaLight-DHc

## Full-range Compact Light Source

### AvaLight-DHc



Get the best out of two worlds with the AvaLight-DHc. It has both a deuterium light source and a halogen light source, providing you with adequate light between 200 and 2500 nm for nearly all absorbance chemistry applications. Deuterium emits light between 200 and 550 nm, where the halogen takes over up to 2500 nm. Coupling this lightsource to the rest of your spectroscopy system is easy with the SMA connector.

This light source is recommended in settings with large fiber cables or direct-attachment to a cuvette holder such as the CUV-DA, due to its relatively low output energy. The integrated TTL-shutter makes saving a dark measurement very simple in combination with AvaSoft (extra IC-DB26-2 needed).

Optionally the AvaLight-DHc is available in a rack-mountable version, to be used in the 19" rack or the 9.5" desktop system.

- \* Combined Deuterium-Halogen
- \* Low cost
- \* Integrated TTL-shutter

A direct-attach cuvette holder CUV-DA (see section accessories) is available for fluorescence or absorbance measurements.




Figure 8 Spectral output of AvaLight-DHc

### Technical Data

	Deuterium Light Source	Halogen Light Source
<b>Wavelength Range</b>	200 - 400 nm	400 - 2500 nm
<b>Stability</b>	< 1 mAU	< 1 mAU
<b>Warm-up time</b>	8 min	1 min
<b>Drift</b>	< 0.25% / h	< 0.25% / h
<b>Optical Power in 600 <math>\mu</math>m fiber</b>	0.2 $\mu$ Watt	7 $\mu$ Watt
<b>Lamp Lifetime</b>	1000 hours	2000 hours
<b>Temp. Range</b>	5°C - 35°C	
<b>Power Supply</b>	12VDC / 450 mA	
<b>Dimensions</b>	175 x 110 x 44 mm	
<b>Lifetime shutter</b>	1.000.000 cycles (typical)	

# AvaLight-HAL-(S)-MINI Tungsten-Halogen Light Source

## AvaLight-HAL-S-MINI



From visible light to near infrared, that's where the AvaLight-HAL-Mini works best. It's a compact, stabilized halogen light source, with adjustable focusing of the fiber connection, maximizing output power at the desired wavelength. The light source also has adjustable output power to provide extra power or longer bulb life.

A filter-slot mounted on the front of the AvaLight-HAL-Mini accepts 1" round or 2" x 2" square filters, to block specific ranges of wavelengths or instantly lower the intensity.

The adjustable focus on the AvaLight-HAL-Mini helps you get the most out of your light source: it makes sure all possible power is transmitted through your optical fiber. Bulb replacement is easy and can be done in a matter of minutes.

Optionally a combined direct-attach cuvette holder and attenuator is available (CUV-ATT-DA-HAL). for attenuation you can use the Inline Filterholder, FH-INL, or the Inline attenuator, ATT-INL. The optical output can be controlled through a dongle at the backside or from your spectrometer. At low setting the

lamp has a color temperature of 2700K but provides over 13000 hours of lifetime. The standard or medium setting changes the color temperature to 2850K and provides 50% more power with a bulb lifetime of 4000 hours. The high power setting gives a color temperature of 3000K, double power compared to the long-life setting and gives you up to 1000 hours of lifetime.

The AvaLight-HAL-S-Mini features an internal TTL-shutter, controllable from your AvaSpec spectrometer. This gives you the ability to use the auto-save dark option in AvaSoft spectroscopy software.




Figure 9 Spectral output of AvaLight-HAL-Mini

## Technical Data

	AvaLight-HAL-(S)-Mini (standard)	AvaLight-HAL-(S)-Mini (long life)	AvaLight-HAL-(S)-Mini (high power)
<b>Wavelength Range</b>	360-2500 nm		
<b>Stability</b>	± 0.1%/ °C		
<b>Time to stabilize</b>	Ca. 10 min.		
<b>Output to bulb</b>	12.0 VDC/ 0.83A	11.3 VDC/ 0.8A	14.1 VDC/ 1.0A
<b>Bulb Life</b>	4000 hrs	> 13000 hrs	< 1000 hrs
<b>Min. Optical power* 200 µm fiber</b>	0.5 mWatt	0.35 mWatt	0.7 mWatt
<b>Min. Optical power* 600 µm fiber</b>	4.5 mWatt	3.2 mWatt	6 mWatt
<b>Min. Optical power* 1000 µm fiber</b>	10 mWatt	7 mWatt	14 mWatt
<b>Bulb Color Temperature</b>	2,850 K	2,730 K	3,000 K
<b>Power requirement</b>	12 VDC / 1.5A		
<b>Temperature range</b>	0-55 °C		
<b>Dimensions, weight</b>	150 x 78 x 37 mm, 510 grams		
<b>Lifetime shutter</b>	1,000,000 cycles (typical)		

\* Optical power measured from 350-1100nm

## Separate 50x50 mm filters to install in AvaLight-HAL-(S)-Mini

<b>GL-WG305-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 305 nm
<b>GL-KG3-3</b>	Separate 50 x 50 x 3 mm band-pass filter, transparent > 325 nm and < 700 nm
<b>GL-BG28-3</b>	Separate 50 x 50 x 3 mm band-pass filter, transparent > 360 nm and < 500 nm
<b>GL-GG385-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 385 nm
<b>GL-GG475-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 475 nm
<b>GL-OG515-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 515 nm
<b>GL-OG550-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 550 nm
<b>GL-OG590-3</b>	Separate 50 x 50 x 3 mm long-pass filter > 590 nm
<b>GL-NG9-1</b>	Separate 50 x 50 x 1 mm Neutral Density filter (transmission 10%, 400-1100 nm)
<b>GL-NG9-2</b>	Separate 50 x 50 x 2 mm Neutral Density filter (transmission 1%, 400-1100 nm)
<b>GL-NG9-3</b>	Separate 50 x 50 x 3 mm Neutral Density filter (transmission 0.1%, 400-1100 nm)

More filter types available, please contact us for ordering information

## Ordering Information

<b>AvaLight-HAL-Mini</b>	• 10W Tungsten Halogen Lamp, fan-cooled, needs extra PS-12V/1.5A power supply
<b>AvaLight-HAL-S-Mini</b>	• As AvaLight-HAL-Mini, incl. TTL shutter, needs extra PS-12V/1.5A power supply
<b>AvaLight-HAL-S-RM</b>	• Rack-mounted version of AvaLight-HAL-S-Mini
<b>IC-DB26-2</b>	• Interface cable AvaSpec-USB2 platform to AvaLight-HAL-(S)-Mini
<b>AvaLight-HAL-B-Mini</b>	• 10W Tungsten Halogen Replacement bulb for AvaLight-HAL-(S)-Mini
<b>PS-12V/1.5A</b>	• Power supply 100-240VAC/12VDC, 1.5A, necessary for AvaLight-HAL-Mini
<b>DONGLE-Mini-H</b>	• Dongle for high power setting
<b>DONGLE-Mini-L</b>	• Dongle for long life setting

# AvaLight-DH-S Deuterium-Halogen Light Source

## AvaLight-DH-S



In need of more power than the AvaLight-DHc? The AvaLight-DH-S is Avantes' most powerful deuterium halogen source. Like the DHc it is also a combined deuterium and halogen light source, capable of transmitting light in the UV/VIS/NIR-range, but has 35 times more halogen output and up to 300 times more deuterium power. The source has a prominent 656 nm deuterium peak which can limit dynamic range (see Avalight-DH-S-BAL as an alternative). It includes a focusing lens assembly, to fully utilize the possibilities and size of your fiber.

The AvaLight-D-S is a deuterium light source only, making it a great option for measurements in the UV range, 190-400 nm. The AvaLight-D-S-DUV version starts even lower at 175 nm, for your deep-UV experiments. This version also offers twice the intensity at 200nm.

The output of the AvaLight-DH-S is optimized for fibers or bundles up to 600 micrometers. For larger fibers the focal point is manually adjustable to optimize the light coupling into your fiber.

The AvaLight-D(H)-S features an integrated TTL-shutter and filter holder for filters of up to 50x50x5.0 mm.

- \* Combined Deuterium-Halogen
- \* UV-VIS-NIR
- \* Deep-UV optional
- \* Powerful




Figure 10 Spectral output AvaLight-DH-S

## Technical Data

	Deuterium (Deep-UV) Long life	Deuterium (Standard) Long life	Halogen
<b>Wavelength Range</b>	175-400 nm	190-400 nm	360-2500 nm
<b>Warm-up Time</b>	30 min.		20 min.
<b>Lamp Power</b>	78W / 0.75A		5W /0.5A
<b>Lamp Lifetime</b>	2000 h		1000 h
<b>Noise (AU)</b>	2x 10 <sup>-5</sup>		10 <sup>-4</sup>
<b>Max. drift</b>	± 0.5%/h		±0.1%/h
<b>Color Temperature</b>	-	-	3000 K
<b>Optical Power* in 200µm fiber</b>	11 µW	11 µW	43 µW
<b>Optical Power* in 600µm fiber</b>	72 µW	72 µW	239 µW
<b>Optical Power* in 1000µm fiber</b>	206 µW	206 µW	354 µW
<b>Power consumption</b>	90 Watt (190Watt for heating D-Lamp 4-5 sec.)		
<b>Power Requirements</b>	100-240VAC 50/60 Hz		
<b>Dimensions / Weight</b>	315 x 165 x 140 mm / ca 5 kg.		
<b>Lifetime shutter</b>	1,000,000 cycles (typical)		

\* total power for the specified wavelength range

For a table of separate 50x50 mm filters to install in AvaLight-D(H)-S see AvaLight-HAL.

## Ordering Information

<b>AvaLight-D-S</b>	• Deuterium light source, 190-400 nm, incl. TTL shutter, -SR fibers needed
<b>AvaLight-DH-S</b>	• Deuterium-Halogen light source, 190-2500 nm, incl. TTL shutter, -SR fibers needed
<b>AvaLight-D-S-DUV</b>	• Deep-UV deuterium light source, 175-400 nm , TTL shutter, -SR fibers needed, 2000h
<b>AvaLight-DH-S-DUV</b>	• Deep-UV deuterium-halogen light source, 175-2500 nm , TTL shutter, -SR fibers needed
<b>IC-DB26-2</b>	• Interface cable AvaSpec-USB2 platform to AvaLight-D(H)S-BAL
<b>AvaLight-D-B</b>	• Replacement deuterium bulb for AvaLight-D/AvaLight DH-BAL light source
<b>AvaLight-D-B-DUV</b>	• Replacement deep-UV deuterium bulb for AvaLight-D/AvaLight DH-BAL light source
<b>AvaLight-DH-B</b>	• Replacement halogen bulb for AvaLight-DH-BAL light source
<b>CUV-DA-DHS</b>	• Direct-attach cuvette holder for AvaLight-D(H)S-BAL

## AvaLight-DH-S-BAL Balanced Power

### AvaLight-DH-S-BAL



The AvaLight-DH-S is a powerful deuterium halogen source, but like any unbalanced deuterium halogen source it does have a very dominant alpha peak at 656 nm. This is why Avantes developed the DH-S-BAL, in which this peak is drastically reduced by a dichroic filter. This means less power, but an increase in the dynamic range of a factor 20. A comparison spectrum, which is taken with a standard AvaSpec-2048, is shown on the next page.

The light source delivers a continuous spectrum with high efficiency. The highest stability is in the ultraviolet, visible and near infrared range, from 215 to 2500 nm. An integrated TTL-shutter and filter holder for filters of up to 50x50x5.0 mm are included. The TTL-shutter can be controlled from any AvaSpec spectrometer, which means the auto-save dark-option in AvaSoft software can be used (please note: IC-DB26-2 cable needed).

Connection to the fiber is done through an SMA-905 connector, which features an adjustable focusing lens assembly. This ensures you getting the maximum possible power into your fiber. For all deuterium light sources solarization resistant fibers (-SR) are recommended (see the fiber-optic section of this catalog). The output of the AvaLight-DH-S-BAL is optimized for fibers or bundles up to 600  $\mu\text{m}$ .

- \* Balanced light source
- \* Wide spectrum: 215-2500 nm
- \* Integrated TTL shutter
- \* High efficiency
- \* Increased dynamic range

The filter holder can be easily replaced by a direct-attach cuvette holder CUV-DA-DHS (see section accessories) useful for fluorescence or absorbance measurements.

	Balanced Deuterium (Standard)	Balanced Halogen Lamp
<b>Wavelength Range</b>	215-500 nm	500-2500 nm
<b>Warm-up Time</b>	30 min.	20 min.
<b>Lamp Power</b>	78 W / 0.75 A	5 W / 0.5 A
<b>Lamp Lifetime</b>	2000 hrs	1000 hrs
<b>Noise (AU)</b>	$2 \times 10^{-5}$	$10^{-4}$
<b>Max. drift</b>	$\pm 0.5\%/hr$	$\pm 0.1\%/hr$
<b>Color Temperature</b>	-	3000 K
<b>Optical Power in 200 <math>\mu\text{m}</math> fiber</b>	6 $\mu\text{W}$	17 $\mu\text{W}$
<b>Optical Power in 600 <math>\mu\text{m}</math> fiber</b>	33 $\mu\text{W}$	160 $\mu\text{W}$
<b>Optical Power in 1000 <math>\mu\text{m}</math> fiber</b>	90 $\mu\text{W}$	448 $\mu\text{W}$
<b>Power consumption</b>	90 Watt (190 Watt for heating D-Lamp 4-5 sec.)	
<b>Power Requirements</b>	100-240VAC 50/60 Hz	
<b>Dimensions / Weight</b>	315 x 165 x 140 mm / ca 5 kg.	
<b>Lifetime shutter</b>	1,000,000 cycles (typical)	




Figure 11 Spectral output AvaLight-DH-S-BAL (red) vs. AvaLight-DH-S (blue)

### Ordering Information

<b>AvaLight-D-S-BAL</b>	• Balanced Deuterium light source, 215-500nm, incl. TTL shutter, -SR fibers needed
<b>AvaLight-DH-S-BAL</b>	• Balanced Deuterium-Halogen light source, 215-2500 nm, incl. TTL shutter, -SR fibers needed
<b>IC-DB26-2</b>	• Interface cable AvaSpec-USB2 platform to AvaLight-D(H)S-BAL
<b>AvaLight-D-B</b>	• Replacement deuterium bulb for AvaLight-D/AvaLight DH-BAL light source
<b>AvaLight-DH-B</b>	• Replacement halogen bulb for AvaLight-DH-BAL light source
<b>CUV-DA-DHS</b>	• Direct-attach cuvette holder for AvaLight-D(H)S-BAL

Add flexibility  
to your spectrometer with  
the Replaceable Slit (-RS) option



# AvaLight-XE Pulsed Xenon

## AvaLight-XE



Perfect for ultraviolet applications like fluorescence, the AvaLight-XE is a pulsed xenon light source. When connected to your AvaSpec spectrometer through the IC-DB26-2 cable (sold separately), the flashes are synchronized with the data collected by the spectrometer. In AvaSoft the number of flashes per scan can be selected.

With a special DUV bulb the AvaLight-XE can be used for deep ultraviolet application (below 200 nm). A special direct-attach cuvette holder is available for your fluorescence applications. For transmission measurements, the AvaLight-XE can be used in conjunction with the CUV-ATT-DA which has an iris attenuator to limit the light output and to avoid saturation.

- \* Pulsed light source
- \* Perfect for fluorescence
- \* Cuvette holder available
- \* Long lifetime




Figure 12 Spectral output of the AvaLight-XE

## Technical Data

<b>Spectral Output</b>	200 nm to 1000 nm
<b>Total Optical Power output</b>	39 $\mu$ J per pulse (average 66 mW)
<b>Optical Power in 200 <math>\mu</math>m fiber</b>	0.66 $\mu$ J per pulse (average 20 $\mu$ W)
<b>Optical Power in 600 <math>\mu</math>m fiber</b>	3.2 $\mu$ J per pulse (average 320 $\mu$ W)
<b>Optical Power in 1000 <math>\mu</math>m fiber</b>	7.4 $\mu$ J per pulse (average 744 $\mu$ W)
<b>Synchronization Input</b>	15 pin sub-D connector, TTL level
<b>Pulse Duration</b>	5 $\mu$ s (at 1/3 height)
<b>Pulse delay</b>	6 $\mu$ s
<b>Pulse rate (max.)</b>	100 Hz
<b>Bulb Life</b>	min. 10 <sup>9</sup> pulses
<b>Connector</b>	SMA-905 connector
<b>Power requirement</b>	12 VDC/550 mA
<b>Dimensions, weight</b>	175 x 110 x 44 mm, 540 grams

## Ordering Information

<b>AvaLight-XE</b>	• Xenon Light Source (200-1000 nm), needs interface cable and power supply
<b>AvaLight-XE-DUV</b>	• DUV (160-1000 nm) version of the AvaLight-XE, needs interface cable and power supply
<b>AvaLight-XE-B</b>	• Spare bulb for the AvaLight-XE (200-1000 nm)
<b>AvaLight-XE-B-DUV</b>	• Spare bulb for the AvaLight-XE-DUV (160-1000 nm)
<b>IC-DB26-2</b>	• Interface cable AvaSpec-USB2 platform to AvaLight-XE
<b>CUV-DA</b>	• Direct-attach cuvette holder for AvaLight-DHc/XE/LED
<b>ATT-DA</b>	• Direct-attach attenuator for AvaLight-DHc/XE/LED
<b>CUV-ATT-DA</b>	• Direct-attach cuvette holder and attenuator for AvaLight-DHc/XE/LED
<b>PS-12V/1.0A</b>	• Power supply 100-240VAC/12VDC, 1.0A for AvaLight-XE