


# MIIPSBox640

**Biophotonic**  
SOLUTIONS INC  
ULTRAFAST PULSE SHAPING & COMPRESSION

## High-Resolution Pulse Shaper

- Measures and compresses pulses in seconds
- Uses the MIIPS® auto-calibration and auto-compression technology
- Complete pulse shaping solution (includes computer, spectrometer and nonlinear optical detection)
- Finalist of the 2008 Prism Award



Push-button interferometric autocorrelation

## Eliminate Manual Tweaking

*With Push-Button Pulse Characterization*

**Includes the MIIPS® technology:**

MIIPS® is an automated procedure for measurement and compression of optical pulses. It uses a calibrated pulse shaper to introduce a set of reference phase functions and monitors their effect on spectrally resolved nonlinear response such as second harmonic generation.

Mathematical analysis of the recorded spectra provides a direct measurement of high-order pulse dispersion. The measured spectral phase can be compensated by the pulse shaper to compress the laser pulses to their transform limit at the target, without manual tweaking.



## System Specifications

SLM type	1D array, liquid crystal
Pixel pitch	100 $\mu\text{m}$
Number of pixels	640
Operating wavelength range	430 - 1700 nm
Maximum spectral window	up to one octave
Maximum average power	2 W
Maximum input pulse energy <sup>&amp;</sup>	250 $\mu\text{J}$
Recommended $1/e^2$ beam diameter	2.5-3.5 mm
Input polarization	linear, horizontal
Dimensions L x W x H	476 x 305 x 260 mm (18.75 x 12.0 x 10.25 in.)
Shaping of spectral phase and amplitude*, independent of the laser repetition rate.	

<sup>&</sup>Depends on the center wavelength of the laser source.

\*Phase-only or phase-and-polarization shaping modes are available.

### Ideal for supercontinuum compression



JOSA B 25, A140 (2008)

Copyright (c) 2010-2016 Biophotonic Solutions, Inc.