

## **KEY FEATURES**

THz.

- COVERS THE ENTIRE THZ SPECTRUM Get the best precision across the entire wavelength range and relative measurements from 30 THz to 0.1
- ROOM TEMPERATURE OPERATION Easier to use and less expensive than a Golay cell.
- > CALIBRATED AT 10.6 μM

THZ-D detectors are calibrated at a single wavelength 10.6 um (30 THz) and at 10 Hz chopping frequency for the THZ9D. Both include typical wavelength correction data from 10.6 to 440 µm. They are used for relative measurements outside that range.

> LARGE AREA

Models range from 9 mm Ø for the THZ9D and 12 mm Ø for the THZ12D.

WIDE RANGE OF MEASUREMENTS

Measure from 100 uW to 3 W of continuous power with the THZ12D model, the highest in our terahertz range of products, and down to 5 uW to 25 mW with the THZ9D model.

USE WITH A UNIVERSAL DISPLAYS & PC INTERFACE

No need for an exclusive monitor. These unique THz detectors work with our display & PC interface.

SDC-500 OPTICAL CHOPPER

The THZ9D model requires the use of an optical chopper, like our SDC-500, running at 10 Hz.

## **OUTPUT OPTIONS**

- SMART DB15 CONNECTOR Contains all the calibration data
- ANALOG OUTPUT When used with APM analog power supply moddule
- > integra ALL-IN-ONE-METER (for THZ12D only) Connects directly to a PC Two models available:
  - USB output (-INT)
  - RS-232 output (-IDR)

## **COMPATIBLE DISPLAYS & PC INTERFACES**



MIRO ALTITUDE



MAESTRO



U-LINK



M-LINK



APM analog power module (for THZ9D)

## **ACCESSORIES**



Stand with steel post (for THZ12D)



Stand with delrin (for THZ9D)



SDC-500 digital optical chopper



Pelican carrying case



Extension cables (4, 15, 20 or 25 m)









	THZ9D-20mS-BL	THZ12D-3S-VP
MAXAVERACE POWER	25mW	3W
EFFECTIVE APERTURE	9 m m 0	12mm0
COMPATIBLE DISPLAYS & PC INTERFACES	MIRO ALTITUDE, MAESTRO, U-LINK, M-LINK &APM	MIRO ALTITUDE, MAESTRO, U-LINK & M-LINK
MEASUREMENT CAPABILITY		
Spectral range•		
Frequency	O.I-30THz	01 -30THz
Wavelength	3000-10 μm	3000-10 µm
Maximum average power		
with MAESTRO	20mW	3W
with U-LINK, M-LINK or MIRO ALTITUDE	25mW	3W
Noise equivalent power b	300nW	O.S <sub>L</sub> W
Minimum measurable power'	N/A	50-100 μW
Thermal drift	N/A	12 μW/C
Rise time (nominal) '	<0.2s	3s
Minimum repetition rate'	IOOOHz	7Hz
Chopping frequency	10 Hz (required)	N/A
Calibration uncertainty 9	±50% at 10.6 μm; ±19% at 10.6 - 440 μm •	±3.0% at 10.6 μm ±8.0% at 10.6 - 300 μm ±15% at 300 - 440 μm
Repeatability	±0.5%	±0.5%
DAMACE THRESHOLDS		
Maximum average power density h	SO mW/cm'	30W/cm'
Maximum energy density	< Q1 J/cm'	<i cm'<="" j="" td=""></i>
PHYSICAL CHARACTERISTICS		
Effective aperture	9 m m 0	12mm0
Absorber	BL (Black Absorber)	VP (Volume Absorber)
Dimensions	38.10 x 26.2 mm	73H x73Wx28D mm (BOD mm with tube)
Weight (head only)	91g	320g
ORDERING INFORMATION		
Compatible stand	STAND-D-233	STAND-D-233
Product page	□(393)□ 2043244	□ 7776 15 □ Ga <b>46</b> 97 45 6

a  $\,$  From 10 to 440  $\,\mu m$ , spectrometer measurement with multiple laser references validation. From 440 to 600  $\mu\text{m}\text{,}$  spectrometer measurement only. From 600 to 3000 µm, relative measurement only.
This spectral range is subject to change.
b. Nominal value. actual value depends on electrical noise in the measurement system.

フォトテクニカ株式会社

- Actual value depends on ambient conditions and the measurement system.
- With anticipation
- e f.
- Maximum output voltage of sensitivity x maximum power.

  Minimum repetition rate for stable average power measurements.

  Including linearity with power.
- g. Including linearity w h At 1064 nm. I W CW.



Specifications are subject to change without notice T 418.657.8003 I info@gentec-eo.com



FAX:048-871-0068,