POWER MANAGEMENT

NEUTRAL DENSITY (ND) FILTERS - UP TO 1 W


We offer various SM1 threaded absorptive ND (Neutral Density) filters that can be fixed directly on the aperture of the Beamage camera via a SM1 to C-mount adaptor. We also offer SM2 threaded filters that can be fixed on the Beamage-4M-FOCUS via a SM2 to T-Mount adaptor. Subsequent filters can be stacked directly on each other. These filters reduce the intensity of all wavelengths without affecting the wavefront of the beam or distorting the image. Sets of 3 filters or 6 filters as well as individual filters are available. An empty SM1 threaded filter holder is also available for those who would like to use their own ND filters with their camera. It holds 25 mm wide filters.

Each 25 mm filter and each holder comes with a SM1 to C-mount adaptor and each 50 mm filter comes with a SM2 to T-Mount adaptor.

MAIN SPECIFICATIONS

| MODEL | ND0.5 TO ND5.0 |
|--------------------------------|--|
| Spectral Range | 400 nm ^a - 1595 nm |
| Filter Diameter | 25 mm Ø |
| Clear Aperture | 22.5 mm Ø (90% of diameter) |
| Dimensional Tolerance | +0.0/-0.25 mm |
| Optical Density Tolerance | ±5% |
| Parallelism | < 10 arcsec |
| Transmitted Wavefront Error | $<\lambda/10$ at 633 nm |
| Surface Flatness | < \(\lambda / 4 \) |
| Surface Quality | 40-20 Scratch-Dig |
| Maximum Power | 1 W |
| Damage Thresholds | 100 W/cm ² or 3 J/cm ² |

^{*} Data specified at 633 nm



OVERVIEW OF THE MODELS

| 25 MM FILTERS | | 50 MM FILTE | RS | EQUIVALENT ATTENUATION | TRANSMITTANCE @ 633NM | SUBSTRATE |
|-------------------------------------|--------|-------------|--------|---------------------------|--------------------------|-----------|
| MODEL | (P/N) | MODEL | (P/N) | | | |
| ND0.5 | 201094 | ND0.5-F0CUS | 203403 | (1/3,16) | ~32% | NG4 |
| ND1.0 | 201045 | ND1.0-F0CUS | 203404 | (1/10) | ~10% | NG4 |
| ND2.0 | 201046 | ND2.0-F0CUS | 203405 | (1/100) | ~1% | NG9 |
| ND3.0 | 201047 | ND3.0-F0CUS | 203406 | (1/1000) | ~0.1% | NG9 |
| ND4.0 | 202600 | ND4.0-FOCUS | 203407 | (1/10 000) | ~0.01% | NG9 |
| ND5.0 | 202601 | ND5.0-FOCUS | 203408 | (1/100 000) | ~0.001% | NG9 |
| "NDSET-6 | | | | | | |
| (Set of all 6 filters)" | 202605 | - | | See Above | See Above | See Above |
| "NDSET-3 | | | | | | |
| (Set of 3 filters (ND1, ND2, ND3))" | 202606 | - | | See Above | See Above | See Above |
| "ND-H | | | | | | |
| (ND filter holder)" | 202607 | - | | | | |

a. For ND4.0 filter, lower limit with other models.

POWER MANAGEMENT


OPTICAL ATTENUATORS - UP TO 500 W



^{*} Beam profiling camera and stand sold separately, adaptor tube for Beam profiling camera is included


| | BA16-60S | BA16K-150S-H5-D0 | BA16K-500F-H9-D0 | |
|------------------------------------|---------------------------------|---------------------------------|---------------------------------|--|
| AXIMUM POWER | 60 W | 150 W | 500 W | |
| FECTIVE APERTURE | 16 mm Ø | 16 mm Ø | 16 mm Ø | |
| OOLING METHOD | Convection | Convection | Fan | |
| EASUREMENT CAPABILITY | | | | |
| Spectral Range | 200 nm - 2100 nm | 200 nm - 2100 nm | 200 nm - 2100 nm | |
| Integrated Power Meter | None | UP19K-15S-H5-D0 | UP19K-110F-H9-D0 | |
| Fan Input Voltage | N/A | N/A | 12 VDC | |
| Average Attenuation | 1700 @ 1064 nm | 1700 @ 1064 nm | 1700 @ 1064 nm | |
| Maximum Beam Diameter | 16 mm Ø | 16 mm Ø | 16 mm Ø | |
| Optical Wedges Material | UV Fused Silica (uncoated) | UV Fused Silica (uncoated) | UV Fused Silica (uncoated) | |
| Sampled Beam Lateral Shift | 21 mm | 21 mm | 21 mm | |
| Sampled Beam Deviation | 90° | 90° | 90° | |
| Residual Beam Deviation | 4° | 4° | 4° | |
| Polarization correction | Yes (pair of orthogonal wedges) | Yes (pair of orthogonal wedges) | Yes (pair of orthogonal wedges) | |
| IYSICAL CHARACTERISTICS | | | | |
| Aperture Diameter | 16 mm Ø | 16 mm Ø | 16 mm Ø | |
| Dimensions | 45H x 47W x 81L mm | 54H x 50W x 86L mm | 54H x 54H x 126L mm | |
| Weight | 0.26 kg | 0.37 kg | 0.46 kg | |
| RDERING INFORMATION | | | | |
| Product Name | BA16-60S | BA16K-150S-H5-D0 | BA16K-500F-H9-D0 | |
| Product Number (without stand) | 203791 | 203792 | 203793 | |
| Add Extension for INTEGRA (USB) | N/A | -INT | -INT | |
| Add Extension for INTEGRA (RS-232) | N/A | -IDR | -IDR | |
| Add Extension for BLU | N/A | -BLU | -BLU | |

POWER MANAGEMENT



PRESERVES POLARIZATION

The BA Series Optical Attenuators use Fresnel reflection on two optical wedges to pick off a small percentage of the input beam. Since the wedges are oriented orthogonally, the S and P polarization states are switched when the reflected beam hits the second wedge, and the difference in reflectance between the two states cancels out. The incoming beam polarization state and irradiance are thus preserved. The wavefront distortion is negligible and the laser output power stability is not affected.





MONITOR POWER AND PROFILE SIMULTANEOUSLY

The BA16K models include a calibrated thermal power detector that also acts as a beam dump. Simply plug the detector into a Gentec-EO Power Meter to measure and display relative power in real time.

To obtain an absolute measurement of power, you will have to determine a correction factor for the BA16K. This can be accomplished in few simple steps. Note, however, that such a correction factor is dependent on the laser polarization and will only be valid if the polarization is stable over time.

In the near future, a power measurement feature will be added to the PC-Beamage software which will then be able to display the laser power density (i.e. W / cm²).



MODULAR CONCEPT

The "Sampled beam" port can be connected to the Beamage via the included adaptor tube. The adaptor tube is also compatible with our ND filters for additional attenuation (recommended for small beams at high power). An ND4 filter is always included with the purchase of a Beamage profiling camera.

The BA Series Optical Attenuators can also be used, stand-alone, as follows:

- OPTICAL PICK-OFF FOR USE WITH OUR ENERGY OR POWER DETECTORS
- ATTENUATOR FOR OUR HIGH SENSITIVITY DETECTORS LIKE M6, PH, ETC.
- POLARIZATION INSENSITIVE BEAM-SPLITTER WITH NO BACK REFLECTIONS

SEE ALSO

ACCESSORIES FOR BEAM DIAGNOSTICS LIST OF REGULAR ACCESSORIES



〒336-0017 埼玉県さいたま市南区南浦和 1-2-17 7EL: 048-871-0067 FAX: 048-871-0068, http://www.phototechnica.co.jp e-mail:voc@phototechnica.co.jp