Q-TUNE

AIR COOLED, DIODE PUMPED, TUNABLE WAVELENGTH Q-SWITCHED LASER

FEATURES

Seamless laser and Optical Parametric Oscillator (OPO) integration

Turnkey performance due water-free pump laser design

Microprocessor controlled operation with self-optimisation, self-calibration capability

Guaranteed > 2 Gshot pump diode lifetime

Hands-free, automated tuning from 210 to 2300 nm

Up to 5 mJ pulse energy in visible range

<6 cm⁻¹ linewidth

Truly variable up to 100 Hz pulse repetition rate: no performance change from single shot to maximum repetition rate

Internal or external triggering modes

Separate output ports for access to pump laser wavelengths

Low power consumption – from 50 to 150 W depending on model

APPLICATIONS

Temporally resolved laser spectroscopy (for example, Light Induced Fluorescence Spectroscopy (LIFS))


Photo-acoustics imaging


Remote sensing

Metrology

Non-linear laser spectroscopy

Flash Photolysis





Q-TUNE is using Optical Parametric Oscillator (OPO) to produce tunable wavelength in 410 – 2300 nm range with linewidth narrower that 6 cm⁻¹. Optional second harmonic generator extends tuning range to 210 – 410 nm with linewidth narrower that 12 cm⁻¹. Combined with shorter than 5 ns pulse duration and up to 100 Hz pulse repetition rate Q-TUNE is perfect coherent light source for temporally resolved spectroscopy, metrology, photoacoustic imaging, remote sensing.

QLI breakthrough water-free laser crystal pumping technology allows to produce high quality laser beam with up to 100 mJ pulse energy. Advanced laser design resulted in very compact, user-friendly turnkey system that requires little maintenance. There is no chillers or bulky power supplies that one needs to fit under the table. All laser electronics is integrated into housing of the Q-TUNE and the only external module is mains adapter that provides 12 or 27 V DC, 50 – 150 W power (depending on model).

Both pump laser and OPO are controlled trough single Ethernet port via build-in web-server. There is no need to install control software – any computer or even cell phone with modern web-browser will be able to control Q-TUNE. API is also provided for integration with user devices

In addition to tunable wavelength output, the Q-TUNE provides two extra ports for access to pump laser beams. Optional extensions available by request:

- Compact spectrometer for monitoring of OPO wavelength and linewidth.
- Fiber coupled OPO output.

SPECIFICATIONS 1)

MODEL		Q-TUNE				
MODEL	Q-TUNE-C10	Q-TUNE-E10	Q-TUNE-C100			
Wavelength range 2)		410 – 2300 nm				
Pulse repetition rate 3)	10 Hz	10 Hz	100 Hz			
Pulse energy 4)	> 1 mJ	> 5 mJ	> 1 mJ			
Linewidth		< 6 cm ⁻¹				
Pulse duration 5)		< 5 ns				
Pulse-to-pulse stability 6)		< 4.5 % RMS				
Power drift 7)		± 3.0 %				
Polarization	horizontal ir	horizontal in signal range, vertical in idler range				
Typical beam diameter 8)	3 mm	5 mm	3 mm			
Typical beam divergence 9)	< 3 mrad	< 5 mrad	< 3 mrad			
Jitter 10)		< 0.5 ns RMS				
PUMP LASER OUTPUTS 11)						
Max pulse energy at 12)						
Fundamental	15 mJ	50 mJ	15 mJ			
2nd harmonic	7 mJ	20 mJ	7 mJ			
3rd harmonic	5 mJ	20 mJ	5 mJ			
DIMENSIONS						

Laser head (W×L×H)	$390 \times 620 \times 135 \text{ mm}^3$		
Power adapter (W×L×H) 13)	$80 \times 200 \times 60 \text{ mm}^3 \text{ (for +27 V DC output)}$		

OPERATING REQUIREMENTS

135

9.07


Cooling requirements	air cooled (water-free)			
Ambient temperature	15 – 25 °C			
Relative humidity	10 – 80 % (non-condensing)			
Mains voltage	90 - 230 V AC, single phase, 47 - 63 Hz 14)			
Average power consumption	<50 W	<100 W	<150 W	

123




- Due to continuous improvements all specifications are subject to change. Unless stated otherwise all specifications are measured at 450 nm and max pulse repetition rate. The parameters marked typical are not specifications. They are indications of typical performance and will vary with each unit we manufacture.
- Extension down to 210 nm is available with optional second harmonic generator.
- ³⁾ Factory-set pulse repetition rate in internal triggering mode. Pulse repetition rate can be divided-down to 1 Hz.
- ⁴⁾ Measured at 450 nm output. See tuning curves for pulse energies at other wavelengths.
- 5) FWHM level at 1064 nm, measured with 350 ps rise time photodiode.
- 6) Measured during 30 seconds operation after warm-up.
- Over 8 hour period after 20 minutes of warm-up, when ambient temperature variation is less than ±2 °C. Power value is calculated every 1 second.
- 8) Beam diameter is measured 20 cm from laser output at the 4 σ level.
- $^{9)}$ Full angle measured at the 4σ level.
- 10) In respect to falling edge of pump diode triggering pulse.
- ¹¹⁾ Laser pulse energy is optimized for OPO pumping and may differ from stand-alone laser specifications.
- Outputs can be configured for simultaneous or non-simultaneous with OPO operation. Values indicated here are for non-simultaneous operation.
- ¹³⁾ Power adapter size depends of model.
- Laser can be powered from appropriate 12 or 27 VDC power source, depending on model. Please inquire for details.





163



Q-TUNE-E10-SH tuning curve





LT-08412, Vilnius, Lithuania

Phone: +370 5 250 3717