

High Energy Femtosecond OPCPA Systems



UltraFlux FF. Custom high pulse energy femtosecond fixed wavelength laser systems delivering up to 1 J pulse energy with pulse duration down to 10 fs.

High Energy UltraFlux laser series delivers up to **30 TW** peak power operating up to 10 Hz.

Originally built for ELI-ALPS (Extreme Light Infrastructure – Attosecond Light Pulse Source) in Hungary, this series is now available for a wide variety of applications.

The master oscillator is a patented (no. EP2827461 and EP2924500) all-in-fiber Yb doped picosecond laser seed source with two fiber outputs. One seeds the OPCPA Front-End and another seeds the Picosecond Pump Laser. Both outputs originate from the same fiber so they are optically synchronized.

This approach eliminates the need for a complex temporal synchronization system typically present in other OPCPA systems.

The Nd:YAG Picosecond Pump Laser system (PPL) is comprised of several sub-systems: diode pumped Regenerative Amplifier, diode pumped Preamplifier, flash lamp pumped Amplifiers, and Second Harmonic Generators which convert

fundamental 1064 nm wavelength to 532 nm. PPL outputs multiple beams at 532 nm. One beam is directed to NOPCPA Front-End subsystem and others are directed to NOCPA amplification stages.

The Front-End NOPCPA (Non-collinear Optical Parametric Chirped Pulse Amplifier) consists of several sub-systems: Picosecond Optical Parametric Amplifier (ps-OPA) amplifying oscillator output pulses, Grating Compressor compressing ps-OPA output pulses, White Light Generator (WLG) broadening the spectrum of ps-OPA output pulses and Femtosecond Non-collinear Optical Parametric Amplifier (fs-NOPA) amplifying WLG output pulses.

The Stretcher sub-system is a Grism (diffraction gratings combined together with prisms) or Offner type pulse stretcher, which stretches output pulse from NOPCPA Front-End and Dazzler (optional Acousto-Optic Programmable Dispersive Filter) for high order phase compensation.

UltraFlux HE SERIES

FEATURES

- ▶ Based on the novel **OPCPA** (Optical Parametric Chirped Pulse Amplification) technology
- ▶ Patented front-end design (patents no. EP2827461 and EP2924500)
- ▶ Up to **1 J** pulse energy at **5 Hz** repetition rate
- ▶ From **Single Shot** to **100 Hz** pulse repetition rate
- ▶ Down to **10 fs** pulse duration
- ▶ Up to **50 mJ** pulse energy at **100 Hz** repetition rate
 - Excellent pulse energy stability: $\leq 1\%$ RMS
 - Excellent long-term average power stability: $\leq 1.5\%$ RMS over **8-hour** period
- ▶ Perfectly synchronized fs and ps output option available
- ▶ Hands free wavelength tuning
- ▶ High contrast pulses without any additional improvement equipment

APPLICATIONS

- ▶ Broadband CARS and SFG
- ▶ Femtosecond pump-probe spectroscopy
- ▶ Nonlinear spectroscopy
- ▶ High harmonic generation
- ▶ Wake field particle acceleration
- ▶ X-ray generation

Multiple stages of NOPCPA (Non-collinear Optical Parametric Chirped Pulse Amplifiers) are used to amplify the stretched pulse from the Stretcher up to 1 J.

Finally, amplified pulses are compressed back down to fs duration in the Pulse Compressor. Bulk glass compressor (combined together with chirped mirror) or traditional diffraction grating compressor can be used depending on pulse duration required and output energy level.

The built-in Output Diagnostics stage ensures reliable, turn-key operation by monitoring critical parameters such as energy, duration, and beam profile.

SPECIFICATIONS

| Model | UltraFlux FT310 | UltraFlux FT10010 | UltraFlux FF50100-F10 | UltraFlux FF8005 |
|--|------------------------------|------------------------------|-------------------------------------|---|
| MAIN SPECIFICATIONS ¹⁾ | | | | |
| Output energy ²⁾ | | | | |
| Signal | 3 mJ | 100 mJ | 50 mJ | 800 mJ |
| SH output ³⁾ | 0.6 mJ | 3.5 mJ ⁴⁾ | 3.5 mJ ⁴⁾ | 3.5 mJ ⁴⁾ |
| TH output ³⁾ | 150 μ J | 1.2 mJ ⁴⁾ | 1.2 mJ ⁴⁾ | 1.2 mJ ⁴⁾ |
| FH output ³⁾ | 30 μ J | 300 μ J ⁴⁾ | 300 μ J ⁴⁾ | 300 μ J ⁴⁾ |
| Pulse repetition rate | 10 Hz | 10 Hz | 100 Hz | 5 Hz |
| Wavelength tuning range | | | | |
| Signal | 750 – 960 nm | 750 – 960 nm | 840 nm | 840 nm |
| SH output ³⁾ | 375 – 480 nm | 375 – 480 nm | 420 nm | 420 nm |
| TH output ³⁾ | 250 – 320 nm | 250 – 320 nm | 280 nm | 280 nm |
| FH output ³⁾ | 210 – 230 nm | 210 – 230 nm | 210 nm | 210 nm |
| Scanning steps | | | | |
| Signal | 5 nm | 5 nm | – | – |
| SH output ³⁾ | 5 nm | 5 nm | – | – |
| TH output ³⁾ | 3 nm | 3 nm | – | – |
| FH output ³⁾ | 1 nm | 1 nm | – | – |
| Pulse duration ^{5) 6)} | 40 \pm 20 fs | 40 \pm 20 fs | \leq 10 fs | 40 \pm 20 fs |
| Pulse energy stability ⁷⁾ | \leq 1.5 % | \leq 1.5 % | \leq 1 % | \leq 1.5 % |
| Long-term power drift ⁸⁾ | \pm 1.5 % | \pm 1.5 % | \pm 1.5 % | \pm 1.5 % |
| Beam spatial profile | Super-Gaussian ⁹⁾ | Super-Gaussian ⁹⁾ | Super-Gaussian ⁹⁾ | Super-Gaussian ⁹⁾ |
| Beam diameter ¹⁰⁾ | \sim 5 mm | \sim 30 mm | \sim 80 mm | \sim 70 mm |
| Beam pointing stability ¹¹⁾ | \leq 30 μ rad | \leq 30 μ rad | \leq 30 μ rad | \leq 30 μ rad |
| Temporal contrast ¹²⁾ | | | | |
| APFC (within \pm 50 ps) | 10 ¹⁰ : 1 | 10 ¹⁰ : 1 | 10 ¹⁰ : 1 | 10 ¹⁰ : 1 |
| Pre-pulse (\leq 50 ps) | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 |
| Post-Pulse (>50 ps) | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 |
| Optical pulse jitter ¹³⁾ | | | | |
| Trig out | \leq 100 ps | \leq 100 ps | \leq 100 ps | \leq 100 ps |
| Pre-Trig out | \leq 50 ps | \leq 50 ps | \leq 50 ps | \leq 50 ps |
| With –PLL option | \leq 2 ps | \leq 2 ps | \leq 2 ps | \leq 2 ps |
| Polarization | Linear | Linear | Linear | Linear |
| PHYSICAL CHARACTERISTICS ¹⁴⁾ | | | | |
| Laser head size (W×L×H mm) | 900 × 1500 × 300 | 1200 × 2000 × 300 | 1200 × 3600 × 500 | 1500 × 2000 × 500, 2 pc. 1200 × 2500 × 500 |
| Power supply size (W×L×H mm) | 553 × 600 × 850 | 553 × 600 × 1200 | 553 × 600 × 1020 553 × 600 × 500 | 553 × 600 × 1800, 2 pc. 553 × 600 × 500 |
| Umbilical length ¹⁵⁾ | 5 m | 5 m | 2.5 m | 5 m |

| Model | UltraFlux FT310 | UltraFlux FT10010 | UltraFlux FF50100-F10 | UltraFlux FF8005 |
|--|--|--|--|--|
| OPERATING REQUIREMENTS ¹⁶⁾ | | | | |
| Electrical power | 200 – 240 V AC, single-phase, 47 – 63 Hz | 200 – 240 V AC, single-phase, 47 – 63 Hz | 208, 380 or 400 V AC, three-phase, 50/60 Hz ¹⁷⁾ | 208, 380 or 400 V AC, three-phase, 50/60 Hz ¹⁷⁾ |
| Power consumption ¹⁸⁾ | ≤ 1 kVA | ≤ 3.5 kVA | ≤ 6 kVA | ≤ 11 kVA |
| Water supply | ≤ 3 l/min, 2 Bar, max 20 °C | ≤ 6 l/min, 2 Bar, max 20 °C | ≤ 10 l/min, 2 Bar, max 20 °C | ≤ 14 l/min, 2 Bar, max 15 °C |
| Operating ambient temperature | 22 ± 2 °C | 22 ± 2 °C | 22 ± 2 °C | 22 ± 2 °C |
| Storage ambient temperature | 15 – 35 °C | 15 – 35 °C | 15 – 35 °C | 15 – 35 °C |
| Relative humidity (non-condensing) | ≤ 80 % | ≤ 80 % | ≤ 80 % | ≤ 80 % |
| Cleanness of the room | ISO Class 7 | ISO Class 7 | ISO Class 7 | ISO Class 7 |

- ¹⁾ Due to continuous improvement, all specifications are subject to change without notice. The parameters marked 'typical' are indications of typical performance and will vary with each unit we manufacture. Presented parameters can be customized to meet customer's requirements.
- ²⁾ Maximum pulse energy specified at 840 nm, SH output at 420 nm, TH output at 280 nm and FH output at 210 nm.
- ³⁾ Harmonic outputs are optional. Specifications valid with respective harmonic module purchased. Outputs are not simultaneous. Maximum harmonic energy depends on OPCPA signal beam profile and pulse duration.
- ⁴⁾ Maximum pump energy for harmonics is limited to 10 mJ @ 840 nm.
- ⁵⁾ Standard pulse duration changes though the wavelength range – shortest pulse duration is achieved ~840 nm spectral range.
- ⁶⁾ Separate 'F10' option can be ordered to reduce pulse duration to ≤ 10 fs. Wavelength tunability not available with 'F10' option.
- ⁷⁾ Under stable environmental conditions, normalized to average pulse energy (RMS, averaged from 60 s).
- ⁸⁾ Measured over 8 hours period after 30 min warm-up when ambient temperature variation is less than ±2 °C.
- ⁹⁾ Super-Gaussian spatial mode of 6-11th order in near field.
- ¹⁰⁾ Beam diameter is measured at signal output at 1/e² level for Gaussian beams and FWHM level for Super-Gaussian beams.
- ¹¹⁾ Beam pointing stability is evaluated as movement of the beam centroid in the focal plane of a focusing element (RMS, averaged from 60 s).

- ¹²⁾ Pulse contrast is only limited by amplified parametric fluorescence (APFC) in the temporal range of ~90 ps which covers OPCPA pump pulse duration. APFC contrast depends on OPCPA saturation level. Our OPCPA systems are ASE-free and pulse contrast value in nanosecond range is limited only by measurement device capabilities (third-order autocorrelator). There are no pre-pulses generated in the system and post-pulses are eliminated by using wedged transmission optics.
- ¹³⁾ Optical pulse jitter with respect to electrical outputs:
 - Trig out > 3.5 V @ 50 Ω
 - Pre-Trig out > 1 V @ 50 Ω
 - PLL option > 1 V @ 50 Ω
- ¹⁴⁾ System sizes are preliminary and depend on customer lab layout and additional options purchased.
- ¹⁵⁾ Longer umbilical with up to 10 m for flash lamp pumped and up to 5 m for diode pumped systems available upon request.
- ¹⁶⁾ The laser and auxiliary units must be settled in such a place void of dust and aerosols. It is advisable to operate the laser in air conditioned room, provided that the laser is placed at a distance from air conditioning outlets. The laser should be positioned on a solid worktable. Access from one side should be ensured.
- ¹⁷⁾ Voltage fluctuations allowed are +10 % / -15 % from nominal value.
- ¹⁸⁾ Required current rating can be calculated by dividing power rating by mains voltage. Power rating is given in apparent power (kVA) for systems with flash lamp power supplies and in real power (kW) for systems without flash lamp power supplies where reactive power is neglectable.



OPTIONS

| Option | Description | Comment |
|-----------|--|---|
| -F10 | Short Pulse option reduces output pulse duration to ≤ 10 fs | Wavelength tunability not available with 'F10' option |
| -CEP | CEP stabilization to ≤ 400 mrad | Passive and active CEP stabilization |
| -DM | 'Deformable Mirror' option for Strehl ration improvement to > 0.9 | |
| -SH/TH/FH | Second, third and fourth harmonic outputs | Typical conversion efficiency from signal is $\sim 35\%$, $\sim 12\%$ and $\sim 3\%$ respectively and depends on beam profile and pulse duration of the system. Harmonic outputs are not simultaneous with signal output |
| -ps out | Additional narrow spectra ps output that is optically synchronized to main system output | Can be simultaneous and non-simultaneous to the main system output. Offers full optical synchronization to fs pulses |
| -AW | Air-Water cooling | No external water required. Heat dissipation equals total power consumption |

PERFORMANCE




Fig 1. Typical UltraFlux FT310 near field beam profile




Fig 2. Typical UltraFlux FT10010 and FF50100-F10 near field beam profile




Fig 3. Typical UltraFlux FF8005 near field beam profile




Fig 4. Typical temporal contrast of UltraFlux FF10010 system

FEMTOSECOND LASERS




Fig 5. Typical long-term power stability of UltraFlux FF5010-F10 system at 840 nm




Fig 6. Typical tuning curve of UltraFlux FT4010 laser system

UltraFlux HE




Fig 7. Typical energies of UltraFlux FT4010 second, third and fourth harmonic outputs




Fig 8. Typical output spectra of UltraFlux FF5010-F10 system




Fig 9. Typical output spectra of UltraFlux FF5010 system at different wavelengths

OUTLINE DRAWINGS




Fig 10. Typical external view of UltraFlux FF5010-F10 system (actual design might vary)




Fig 11. Typical UltraFlux FF5010-F10 laser system external dimensions

POWER SUPPLY

| Cabinet | Usable height | Height H, mm | Width W, mm | Depth D, mm |
|---------|---------------|----------------------------|-------------|-------------|
| MR-9 | 9 U | 455.5 (519 ¹⁾) | 553 | 600 |
| MR-12 | 12 U | 589 (653 ¹⁾) | 553 | 600 |
| MR-16 | 16 U | 768 (832 ¹⁾) | 553 | 600 |
| MR-20 | 20 U | 889 (952 ¹⁾) | 553 | 600 |
| MR-25 | 25 U | 1167 (1231 ¹⁾) | 553 | 600 |

¹⁾ Full height with wheels.




Fig 12. Typical UltraFlux laser system power supply dimensions (MR rack used depends on the laser model)

ORDERING INFORMATION

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer than 1 hour then laser (system) needs warm up for a few hours before switching on.

UltraFlux (1) (2)(3)-(4)

Model

Any additional options:
See 'Options' table

Fixed or tunable wavelength:
FF → fixed wavelength
FT → tunable wavelength

Pulse repetition rate:
5 → 5 kHz
10 → 10 Hz
100 → 100 Hz

Energy level:
3 → 3 mJ
50 → 50 mJ
100 → 100 mJ
800 → 800 mJ

High repetition rate Tunable Wavelength Femtosecond OPCPA Systems



UltraFlux HR series is a compact high repetition rate tunable wavelength femtosecond laser system which incorporates the advantages of dual output ultrafast fiber laser, solid-state and parametric chirped pulse amplification technologies.

A novel OPCPA front-end technology uses a dual output picosecond fiber laser for seeding both picosecond DPSS pump laser and femtosecond parametric amplifier with a spectrally broadened output.

This approach greatly simplifies the system – excludes femtosecond regenerative amplifier and eliminates the need of pump and seed pulse synchronization while ensuring practically zero jitter between the channels. In addition to that, contrast of the output pulses in picosecond to nanosecond time scale is enhanced.

All UltraFlux series laser systems are assembled on a rigid breadboard or optical table to ensure excellent long-term stability. Modular internal design offers high level of customization and easy scalability. All of these systems can be customized according to customer requirements by adding custom specifications or multiple channels.

Incorporation of parametric chirped pulse amplification technology together with a novel ultrafast fiber laser helped to create and bring to the market a new tool for femtosecond pump-probe, nonlinear spectroscopy, emerging high harmonic generation experiments and other femtosecond and nonlinear spectroscopy applications. With this laser ultrafast science breakthrough is closer to any photonics lab than ever before.

UltraFlux HR SERIES

FEATURES

- ▶ Based on the novel OPCPA (Optical Parametric Chirped Pulse Amplification) technology
- ▶ Patented front-end design (patents no. EP2827461 and EP2924500)
- ▶ 750 – 960 nm, 375 – 480 nm, 250 – 320 nm and 210 – 230 nm wavelength tuning ranges
- ▶ Up to 14 mJ pulse energy at 1 kHz repetition rate
 - Excellent pulse energy stability: $\leq 1\%$ RMS
 - Excellent long-term average power stability: $\leq 1.5\%$ RMS over 8-hour period
- ▶ Perfectly synchronized fs and ps output option available
- ▶ Hands free wavelength tuning
- ▶ High contrast pulses without any additional improvement equipment

APPLICATIONS

- ▶ Broadband CARS and SFG
- ▶ Femtosecond pump-probe spectroscopy
- ▶ Nonlinear spectroscopy
- ▶ High harmonic generation

SPECIFICATIONS

| Model | UltraFlux FT031k | UltraFlux FT31k | UltraFlux FT61k | UltraFlux FT141k |
|--|--|--|--|--|
| MAIN SPECIFICATIONS ¹⁾ | | | | |
| Output energy ²⁾ | | | | |
| Signal | 300 µJ | 3 mJ | 6 mJ | 14 mJ |
| SH output ³⁾ | 60 µJ | 0.6 mJ | 1.5 mJ | 3.5 mJ ⁴⁾ |
| TH output ³⁾ | 15 µJ | 150 µJ | 0.4 mJ | 1.2 mJ ⁴⁾ |
| FH output ³⁾ | 3 µJ | 30 µJ | 100 µJ | 300 µJ ⁴⁾ |
| Pulse repetition rate | 1 kHz | 1 kHz | 1 kHz | 1 kHz |
| Wavelength tuning range | | | | |
| Signal ⁵⁾ | 750 – 960 nm | 750 – 960 nm | 750 – 960 nm | 750 – 960 nm |
| SH output ³⁾ | 375 – 480 nm | 375 – 480 nm | 375 – 480 nm | 375 – 480 nm |
| TH output ³⁾ | 250 – 320 nm | 250 – 320 nm | 250 – 320 nm | 250 – 320 nm |
| FH output ³⁾ | 210 – 230 nm | 210 – 230 nm | 210 – 230 nm | 210 – 230 nm |
| Scanning steps | | | | |
| Signal | 5 nm | 5 nm | 5 nm | 5 nm |
| SH output ³⁾ | 5 nm | 5 nm | 5 nm | 5 nm |
| TH output ³⁾ | 3 nm | 3 nm | 3 nm | 3 nm |
| FH output ³⁾ | 2 nm | 2 nm | 2 nm | 2 nm |
| Pulse duration ^{5) 7)} | 40 ± 20 fs | 40 ± 20 fs | 40 ± 20 fs | 40 ± 20 fs |
| Pulse energy stability ⁸⁾ | ≤ 1.5 % | ≤ 1 % | ≤ 1 % | ≤ 1 % |
| Long-term power drift ⁹⁾ | ± 1.5 % | ± 1.5 % | ± 1.5 % | ± 1.5 % |
| Beam spatial profile | Gaussian | Super-Gaussian ¹⁰⁾ | Super-Gaussian ¹⁰⁾ | Super-Gaussian ¹⁰⁾ |
| Beam diameter ¹¹⁾ | ~ 2 mm | ~ 5 mm | ~ 7 mm | ~ 15 mm |
| Beam pointing stability ¹²⁾ | ≤ 30 µrad | ≤ 30 µrad | ≤ 30 µrad | ≤ 30 µrad |
| Temporal contrast ¹³⁾ | | | | |
| APFC (within ± 50 ps) | 10 ¹⁰ : 1 | 10 ¹⁰ : 1 | 10 ¹⁰ : 1 | 10 ¹⁰ : 1 |
| Pre-pulse (≤ 50 ps) | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 |
| Post-Pulse (> 50 ps) | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 | 10 ⁸ : 1 |
| Optical pulse jitter ¹⁴⁾ | | | | |
| Trig out | ≤ 100 ps | ≤ 100 ps | ≤ 100 ps | ≤ 100 ps |
| Pre-Trig out | ≤ 50 ps | ≤ 50 ps | ≤ 50 ps | ≤ 50 ps |
| With -PLL option | ≤ 2 ps | ≤ 2 ps | ≤ 2 ps | ≤ 2 ps |
| Polarization | Linear, Horizontal | Linear, Horizontal | Linear, Horizontal | Linear, Horizontal |
| PHYSICAL CHARACTERISTICS ¹⁵⁾ | | | | |
| Laser head size (W×L×H mm) | 750 × 1200 × 300 | 900 × 1500 × 300 | 900 × 1800 × 300 | 1200 × 2000 × 300 |
| Power supply size (W×L×H mm) | 553 × 600 × 850 | 553 × 600 × 850 | 553 × 600 × 850 | 553 × 600 × 1250 |
| Umbilical length ¹⁶⁾ | 2.5 m | 2.5 m | 2.5 m | 2.5 m |
| OPERATING REQUIREMENTS ¹⁷⁾ | | | | |
| Electrical power | 200 – 240 V AC, single-phase, 47 – 63 Hz | 200 – 240 V AC, single-phase, 47 – 63 Hz | 208, 380 or 400 V AC, three-phase, 50/60 Hz ¹⁸⁾ | 208, 380 or 400 V AC, three-phase, 50/60 Hz ¹⁸⁾ |
| Power consumption ¹⁹⁾ | ≤ 1 kW | ≤ 2 kW | ≤ 5 kW | ≤ 8 kW |
| Water supply | not required | not required | not required | ≤ 5 l/min, 2 Bar, max 20 °C |
| Operating ambient temperature | 22 ± 2 °C | 22 ± 2 °C | 22 ± 2 °C | 22 ± 2 °C |
| Storage ambient temperature | 15 – 35 °C | 15 – 35 °C | 15 – 35 °C | 15 – 35 °C |
| Relative humidity (non-condensing) | ≤ 80 % | ≤ 80 % | ≤ 80 % | ≤ 80 % |
| Cleanness of the room | ISO Class 7 | ISO Class 7 | ISO Class 7 | ISO Class 7 |

¹⁾ Due to continuous improvement, all specifications are subject to change without notice. The parameters marked 'typical' are indications of typical performance and will vary with each unit we manufacture. Presented parameters can be customized to meet customer's requirements.

²⁾ Maximum pulse energy specified at 840 nm, SH output at 420 nm, TH output at 280 nm and FH output at 210 nm.

³⁾ Harmonic outputs are optional. Specifications valid with respective harmonic module purchased. Outputs are not simultaneous. Maximum harmonic energy depends on OPCPA signal beam profile and pulse duration.



- 4) Maximum pump energy for harmonics limited to 10 mJ @ 840 nm.
- 5) Optional extended tuning range of 700 – 1010 nm available upon request.
- 6) Standard pulse duration changes though the wavelength range – shortest pulse duration is achieved ~840 nm spectral range.
- 7) Separate 'F10' option can be ordered to reduce pulse duration to ≤ 10 fs. Wavelength tunability not available with 'F10' option.
- 8) Under stable environmental conditions, normalized to average pulse energy (RMS, averaged from 60 s).
- 9) Measured over 8 hours period after 30 min warm-up when ambient temperature variation is less than ± 2 °C.
- 10) Super-Gaussian spatial mode of 6-11th order in near field.
- 11) Beam diameter is measured at signal output at $1/e^2$ level for Gaussian beams and FWHM level for Super-Gaussian beams.
- 12) Beam pointing stability is evaluated as movement of the beam centroid in the focal plane of a focusing element (RMS, averaged from 30 s).
- 13) Pulse contrast is only limited by amplified parametric fluorescence (APFC) in the temporal range of ~90 ps which covers OPCPA pump pulse duration and is better than $10^7:1$. APFC contrast depends on OPCPA saturation level. Our OPCPA systems are ASE-free and pulse contrast value in nanosecond range is limited only by measurement device capabilities (third-order autocorrelator). There are no pre-pulses generated in the system and post-pulses are eliminated by using wedged transmission optics.
- 14) Optical pulse jitter with respect to electrical outputs:
 - Trig out > 3.5 V @ 50 Ω
 - Pre-Trig out > 1 V @ 50 Ω
 - PLL option > 1 V @ 50 Ω
- 15) System sizes are preliminary and depend on customer lab layout and additional options purchased.
- 16) Longer umbilical with up to 10 m for flash lamp pumped and up to 5 m for diode pumped systems available upon request.
- 17) The laser and auxiliary units must be settled in such a place void of dust and aerosols. It is advisable to operate the laser in air conditioned room, provided that the laser is placed at a distance from air conditioning outlets. The laser should be positioned on a solid worktable. Access from one side should be ensured.
- 18) Voltage fluctuations allowed are +10 % / -15 % from nominal value.
- 19) Required current rating can be calculated by dividing power rating by mains voltage. Power rating is given in apparent power (kVA) for systems with flash lamp power supplies and in real power (kW) for systems without flash lamp power supplies where reactive power is neglectable.

OPTIONS

| Option | Description | Comment |
|-----------|---|---|
| -F10 | Short Pulse option reduces output pulse duration to ≤ 10 fs | Wavelength tunability not available with 'F10' option |
| -CEP | CEP stabilization to ≤ 400 mrad | Passive and active CEP stabilization |
| -DM | 'Deformable Mirror' option for Strehl ration improvement to > 0.9 | |
| -SH/TH/FH | Second, third and fourth harmonic outputs | Conversion efficiency from signal respectively ~20 %, ~5 % and ~1 %. Harmonic outputs are not simultaneous with signal output |
| -ps out | Additional ps output that is optically synchronized to main system output | Can be simultaneous and non-simultaneous to the main system output |
| -AW | Air-Water cooling | No external water required. Heat dissipation equals total power consumption |

PERFORMANCE




Fig 1. Typical UltraFlux FT031k near field beam profile




Fig 2. Typical UltraFlux FT31k near field beam profile




Fig 3. Typical UltraFlux FT61k and FT141k near field beam profile




Fig 4. Long-term power stability measurement at 800 nm wavelength




Fig 5. Typical temporal contrast of UltraFlux systems




Fig 6. Typical energy tuning curve of UltraFlux FT31k laser system




Fig 7. Typical energy tuning curves of UltraFlux FT31k second, third and fourth harmonic outputs




Fig 8. Typical output spectra of UltraFlux FT31k system at multiple wavelengths



Fig 9. Typical external view of UltraFlux FT031k system (actual design might vary)

OUTLINE DRAWINGS




Fig 10. Typical UltraFlux FT031k laser system external dimensions

POWER SUPPLY

| Cabinet | Usable height | Height H, mm | Width W, mm | Depth D, mm |
|---------|---------------|----------------------------|-------------|-------------|
| MR-9 | 9 U | 455.5 (519 ¹⁾) | 553 | 600 |
| MR-12 | 12 U | 589 (653 ¹⁾) | 553 | 600 |
| MR-16 | 16 U | 768 (832 ¹⁾) | 553 | 600 |
| MR-20 | 20 U | 889 (952 ¹⁾) | 553 | 600 |
| MR-25 | 25 U | 1167 (1231 ¹⁾) | 553 | 600 |

¹⁾ Full height with wheels.




Fig 11. Typical UltraFlux laser system power supply dimensions (MR rack used depends on the laser model)

ORDERING INFORMATION

Note: Laser must be connected to the mains electricity all the time. If there will be no mains electricity for longer than 1 hour then laser (system) needs warm up for a few hours before switching on.

UltraFlux (1) (2)(3)-(4)

Model

Any additional options:
See 'Options' table

Fixed or tunable wavelength:
FF → fixed wavelength
FT → tunable wavelength

Pulse repetition rate:
1k → 1 kHz

Energy level:
03 → 300 μJ
3 → 3 mJ
6 → 6 mJ
14 → 14 mJ

PHOTO TECHNICA www.phototechnica.co.jp
フォトテクニカ株式会社
 〒336-0017 埼玉県さいたま市南区南浦和 1-2-17
 TEL:048-871-0067 FAX:048-871-0068
 e-mail:voc@phototechnica.co.jp