

BEAMAGE-M2


Automated M² measurement system



KEY FEATURES

- **LARGE APERTURES**
The only M² system on the market equipped with a complete set of 50mm (2") optics. Also, the sensor is 11.3 x 11.3mm
- **SIMPLE ALIGNMENT**
Two beam-steering mirrors are included for quick and easy alignment of your laser into the system.
- **COMPACT**
The low-profile ingenious mechanics make it easy to fit the device on any optical table
- **ISO COMPLIANT**
The calculations are fully compliant to the ISO 11146 and 13694 standards
- **FAST ACQUISITION**
Make a complete, ISO-compliant measurement in only 20 seconds with the ROI feature and in less than a minute with full-frame acquisition
- **FLEXIBLE & INTUITIVE SOFTWARE**

AUTOMATED MEASUREMENTS



Inside the BEAMAGE-M2, a computer-controlled motorized rail allows precise positioning of two mirrors, which in turn allow a 400 mm beam path difference. At each position of the translation stage, a beam profile is acquired and the beam diameter is measured. The automation of the translation stage allowed by the software is the key to a fast measurement.

PRACTICAL ALIGNMENT TOOL



Each BEAMAGE-M2 system includes an alignment tube that helps you set up the system faster. Simply use the two alignment mirrors to center your laser beam onto both irises, and you will be ready to start measuring in no time!

The fluorescent material around the pinholes also helps to align beams that are in the NIR range without having to use an IR viewer.

BEAMAGE-M2

Specifications



| BEAMAGE-M2 | |
|---|---|
| SENSOR TECHNOLOGY | Beamage-4M included |
| EFFECTIVE APERTURE | Ø 48 mm optics - 11.3 x 11.3 mm sensor |
| MEASUREMENT CAPABILITY | |
| System wavelength range | 350 - 1100 nm |
| Attenuation range | 3 Flip-mount attenuators for 8 levels of attenuation: no attenuation, ND0.5, ND1, ND2, ND1.5, ND2.5, ND3, ND3.5 |
| Beam diameter range ^a | 55 µm to 11.3/3 mm |
| Translation stage | |
| Mechanical travel range | 200 mm |
| Effective optical path range | 400 mm |
| Lens focal length | 5 AR-coated lenses included: 200 mm, 250 mm, 300 mm, 400 mm and 500 mm |
| Typical M ² accuracy ^b | ± 5% |
| Typical M ² repeatability ^b | ± 2% |
| Applicable light sources | CW and quasi-CW |
| Typical measurement time | 45 s with full-frame acquisition |
| DAMAGE THRESHOLDS ^c | |
| Maximum average power | 1 W with ND filter |
| Maximum density (1064 nm) | CW: 10 W/cm ² ; Pulsed: 0.1 J/cm ² |
| PHYSICAL CHARACTERISTICS | |
| Dimensions | |
| Main enclosure | 357 mm (L) x 165 mm (W) x 135 mm (H) |
| Total (including external mirrors) | 602 mm (L) x 193 mm (W) x 172 mm (H) |
| Optical axis height | 86 mm |
| Weight | 6.6 kg |
| Power supply | 48 VDC, 1.25A out |
| SOFTWARE | |
| Displays | 2D, 3D, XY, Beam Tracking and M ² |
| Beam diameter definitions | D4σ 1/e ² along crosshairs (13.5%) FWHM along crosshairs (50%) Custom (%) |
| Beam quality definitions | Laser beam quality M ² : M ² _x , M ² _y (ISO compliant) Beam Propagation Factor: BPP _x , BPP _y Width at waist: W _x , W _y Waist location and offset: Z _x , Z _y , ΔZ Divergence angle: θ _x , θ _y Rayleigh length: Z _{Rx} , Z _{Ry} Astigmatism |
| Printing and reports | Full report in print-ready format |
| ORDERING INFORMATION | |
| Product page |  |

Specifications in the table above are for the use with a Beamage-4M beam profiler (included in the Beamage-M2 kit)

- a. At the Beamage sensor
- b. Depending on the beam quality and optical configuration
- c. With ND4 filter at the Beamage


Specifications are subject to change without notice
T 418.651.8003 | info@gentec-eo.com

BEAMAGE-M2

Software features

Select which set of Rayleigh range boundaries to display on the graph: X, Y or both

M2 Results tab:
View and understand all the measured parameters quickly, for both the initial laser beam and the beam inside the BEAMAGE-M2 system



Switch to 3D or 2D displays to see each of the measured profiles

M2 Setup tab:
Control your acquisition parameters

BEAMAGE-M2

Software features




ISO COMPLIANCE MADE SIMPLE

With the "RUN M2 SETUP" button, the software automatically defines new parameters for a more precise M^2 measurement. The "ISO SCAN" data set complies with the ISO-11146 M^2 measurement standard, being spread between $-3Z_R$ and $+3Z_R$.

The automatic settings are updated after each calculation, considering the values of Z_0 and Z_R from the latest measurement.

By default, the results graph always shows the calculated positions of the first three Rayleigh distances on each side of the waist.



FULL CONTROL ON YOUR DATA

During an M^2 scan, each of the measured profiles is saved and the flexible software gives you complete control on your acquired data.

- View each acquired profile in 2D display or 3D display.
- Add measurement points to a data set at the position of your choice with the "ADD" button.
- Remove unwanted profiles from your data set & recalculate the measurements.



FAST ATTENUATION

Add or remove attenuation with the flick of a finger. The software adjusts the exposure time at each frame during an acquisition, and advises the user on the required attenuation.