Raman Spectroscopy Software

Included with the AvaRaman systems, AvaSoft-Raman enables full control over your Raman spectroscopy system.

In addition to most of the features available in AvaSoft-FULL, AvaSoft-Raman, which is a standalone application, also features:


- Display of the wavelength axis in cm⁻¹
- Auto calibration routines to determine the excitation laser peak (please note that an AvaRaman-Calibration tile is needed, sold separately)
- Integration time progress bar to indicate integration time status for longer spectral acquisitions
- View signal in normalized counts
- Software baseline correction for fluorescence suppression

AvaSoft-Raman also features history channel functions to monitor peak values or concentration versus time.

Process control and export to Excel add on modules are also available for online analyses and control.



Image: Specify and Specify and



Ordering Information

Avasoft-Raman

• Raman application add-on software, to be ordered with AvaSoft-Full

Avasoft-All • Full version AvaSpec software (includes AvaSoft-COL, -IRRAD, -Raman and, -ThinFilm)

Three years limited warranty on all Avantes spectrometers, light sources and accessories



www.phototechnica.co.jp フォトテクニカ株式会社

〒336-0017 埼玉県さいたま市南区南浦和 1-2-17 TEL:048-871-0067 FAX:048-871-0068 e-mail:voc@phototechnica.co.jp

info@avantes.com | www.avantes.com | 79